

# revolve<sup>TM</sup>

*Join the Revolution.*



# Why join the Revolution?

Introducing Revolve™ handpieces by ProDentUSA™ - the ultimate machine for practitioners seeking pinnacle performance and extraordinary results. Crafted with precision and designed for durability, these handpieces are the peak of cutting-edge technology. All Revolve™ handpieces are designed to last and guaranteed to satisfy.

Discover our Revolve™ collection:

## High Speed Handpieces

- Fiber optic
- Lightweight, Titanium casing ensures durability
- Eternal grip, with non-slip coating
- While stainless steel ball bearings expand under high temperatures, wearing out the handpiece and decreasing speeds, Revolve™'s ceramic ball bearings reduce wear and keep temperatures low and don't expand with heat, resulting in increased speed and power.
- Most powerful chuck on the market
- 4-port spray for uniform water output
- Up to 420,000rpm
- 25 Watt power
- Weighs only 1.8oz
- 1-year warranty on turbine, 2-year warranty on casing

## Slow Speed Handpieces

- Up to 20,000rpm
- Reversible motor
- Quiet and smooth motor
- Weighs 2.8oz
- 5-year warranty

Revolve™ handpieces are designed to function flawlessly at your practice. Experience guaranteed, unwavering power when you need it most. Each handpiece is meticulously crafted and diligently tested prior to shipping. We invite you to explore our catalog of the revolutionary handpiece.



Pictured here:  
(Above) Torque-head high speed  
handpiece, R750TK-L

(Below) Mini-head high speed  
handpiece, R550TK-L



Shop  
Revolve™  
online!



## High Speed Handpieces

Revolve™ high speed handpieces come with a 2-year warranty on the casing and 1-year warranty on the turbine. All handpieces are available in either NSK- or Kavo-type backing.



R750GQ  
Self-generating LED,  
4-hole **\$310**  
Non-Optic - R750NQ **\$290**



R750NQL  
Lubricant-free  
Torque head **\$450**



550TK-L  
Mini head size,  
Kavo-style **\$500**  
(550TN-L for NSK-style)



750TK-L  
Torque head size,  
Kavo-style **\$500**  
(750TN-L for NSK-style)

## Surgical Handpieces

Revolve™ surgical handpieces are designed to push air away from the surgical site. The convenient 45° angle is optimal for increased visibility during surgical procedures.





## Electric Handpieces

Revolve™ electric handpieces come with a 1-year warranty and are compatible with any electric motor.



RAD15  
1:5 Increasing Contra-Angle,  
Non-Optic Handpiece \$490



RAD15L  
1:5 Increasing Contra-Angle,  
Fiber-Optic Handpiece \$570



RAD20  
20:1 Implant Angle,  
Non-Optic Handpiece \$620



RAD20L  
20:1 Implant Angle,  
Fiber-Optic Handpiece \$670

## Slow Speed Handpieces



RCALS  
Contra-angle  
Latch set  
\$220

R304-NCAA  
Slow-speed set:  
Slow-speed motor, Nose  
Cone, & Contra-angle  
\$600



R304  
Slow-speed,  
motor only  
\$475

## Contra-Angle Accessories



RLH  
Latch-button head  
\$120



RPB  
Push-button head  
\$175



R310  
Contra-angle  
sheath  
\$100



RCAPBS  
Contra-angle  
Push-button set  
\$270



## Couplers

Revolve™ couplers are compatible with all handpieces and come in Kavo or NSK-styles with 4- or 6-hole tubing configurations.



RLN06  
NSK-Style 6-Hole  
\$150



R4NSKG  
NSK-Style 4-Hole Coupler,  
Self-generating LED  
\$150



RAD-K4  
Kavo-Style 4-Hole  
\$150



RK6H  
Kavo-Style 6-Hole  
\$150



## Join the Revolution.



"These handpieces are better than Kavo, better than NSK! They're insanely light - unparalleled torque. You've gotta get the word out about these!"


- Dr. Cary Wagner, NYC

"Finally! After having searched high and low for an alternative to the Midwest, NSK, Bien Air monopoly I took a chance on these guys. I used to get the AZdent and other super econ handpieces but they went bad within a few months and otherwise would have the top lug nut come loose. These feel and operate like my \$\$\$ handpieces. Tons of torque, half the cost and really solid. I love them."

- Dr. David Donelson, TN

"They are amazing! The price is reasonable and the quality is excellent. The self-generating LED light is very bright and reliable and the turbine runs smoothly. They are well built."

- Dr. James Ford, OH



Ceramic ball bearings prevent power loss and wear by keeping temperature cool



Advanced 4-port spray



Most powerful chuck on the market



27 Watt power



Lightweight, Titanium body



1-year warranty on turbine,  
2-year warranty on casing

